USER MANUAL
Shower and commode chair
Product code: Movano

1. Characteristics
The shower and commode chair is made of durable PP (polypropylene) material. It has height-adjustable aluminum legs finished with non-slip rubber tips. The adjustment allows for customization to meet the individual needs of the patient. The backrest and seat are equipped with a soft, removable cushion. Side armrests increase stability and safety of use. The chair is equipped with an easily extendable plastic sanitary container with a lid.

Set includes: Seat frame, backrest, armrests (2 pcs), legs with tips (4 pcs), seat cushion, backrest cushion, sanitary container with lid.

Elements of the commode chair:

- **A.** backrest
- **B.** backrest cushion
- **C.** armrest
- **D.** seat cushion
- **E.** seat frame
- **F.** leg
- **G.** anti-sliding tip
- **H.** sanitary container

2. Application
The shower and commode chair is an auxiliary equipment that facilitates the user to perform hygienic activities. It is recommended for people who have problems with independent movement, with dysfunctions and movement limitations, which prevents or hinders independent performance of activities of

daily life. The chair has a toilet function thanks to the provision of a sanitary receptacle under the seat. It facilitates carrying out the user's daily physiological and care activities. The chair is also recommended for use in showers and baths, in home and hospital settings.

3. Assembly, adjustment and usage

UNFOLDING

Step 1: Turn the chair frame upside down and insert the legs into the mounting holes. The legs have locking pins - push the pin in before inserting the leg, then lock the leg.

Step 2. Adjust the height of the legs with the adjustment pins.

Step 3: Turn the chair upright and place it on the ground.
Note: All four legs must be adjusted to the same height and rest firmly on the ground for stability.

Step 4: Insert the backrest into the mounting holes in the seat frame (mounting latch).

Step 5: Insert both armrests into the mounting holes in the seat frame (mounting latch).

Step 6: Insert the sanitary container with lid into the hole in the seat frame.

Step 7: Insert the backrest cushion into the mounting holes in the backrest (3 mounting clips).

Step 8: Place the seat cushion on the seat of the chair - the shape of the cushion makes it easy to fit.

HEIGHT ADJUSTMENT

To adjust the height of the chair, press the leg adjustment pins, select the desired height and then lock the set height with the pin. The legs are adjustable in increments of 2.5 cm.
Note: all legs must be set at the same level and stand firmly on the ground.

USAGE

Lift the seat cushion, and then remove the lid of the sanitary container. After use, cover the container with the lid and place the seat cushion on the seat.
Note! Remember to empty the toilet sanitary container after each use.

4. Technical data

Overall width	56 cm
Overall depth	48 cm
Backrest height	35 cm
Seat height from the ground	41-53 cm
Overall height from the ground	79-88 cm
Seat dimensions	46 cm x 43 cm
Weight	5,1 kg
Maximum capacity	150 kg
Colour	grey

5. Notes

Timago International Group does not bear responsibility of improper use of the product, failure to comply safety regulations and misuse.

6. Maintenance

Carry out regular maintenance of the technical condition of the chair:

- Clean the commode chair regularly using only water and mild detergents. Dry thoroughly.
- **Do not use** for cleaning: strong polishing agents and detergents, abrasive agents bleach, petrol or solvents (it may cause frame damage, visible defects and in consequence - loss of warranty rights).
- Empty and disinfect the sanitary container and the toilet seat after each use.

7. Storage and transport conditions

The product should not be exposed to direct sunlight, water or moisture.

8. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

9. Labels

REF	Reference number		Manufacturer
LOT	LOT number	UDI	UDI code
SN	Serial number	MD	Medical device
	Note		Please read the instruction
CE	The product meets the current EU legal requirements for medical devices (MDR 2017/745)		

Our complete line of products is available from our distributors. Learn more at: **www.timago.com**. Thank you for choosing Timago!

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